



Molift Partner 205



Molift Partner 205 handles all moving and handling situations, at home and in institutional environments

- **excellent lifting height**
- **lifts from the floor without height adjustment**
- **unique quick-release system makes it easy and quick to change between various suspensions**

molift[®]
designed for life

Molift Partner 205 makes life easier

- **user friendly**
 - excellent lifting height
 - lifts from the floor without height adjustment
 - low weight and optimum wheel positioning makes Molift Partner 205 easy to manoeuvre
 - 4-point suspension for open seating position and with lots of space around the head
 - the new Molift Stretcher can be mounted directly on the standard Molift Partner 205 4-point suspension
 - **Accessories:** Molift Scale, Molift Stretcher and supporting Arms (optional)
- **4-point suspension**
 - 4-point suspension prevents the user from feeling squeezed, but gives instead an open sitting position with lots of space around the head.
 - 4-point suspension together with a Molift sling gives the user correct and comfortable position from Sitting to lying and vice versa as the loops of the sling slides into a natural position during the lift.
- **automatic service warning feature**
 - Molift Partner 205 is equipped with all necessary safety devices such as emergency stop, and electronic and manual emergency lowering, etc.
 - Electronics register the number of lifts and informs when a service is necessary.



Specification

Lifting capacity	lbs	452
Weight	lbs	83,8
Measures (in)	Lifting stroke	55
	(L x W x H)	48 x 28 x 55
	Leg height	4 ½
Battery	26.4 V NiCd 1.9 Ah / MiMh 2,1 Ah	
Wheels	Back Ø100 mm, lockable, rubber	2
	Front, Ø75 mm, polyurethane	2



5 year warranty when subject to an annual service contract. 1 year on batteries and hand control.

Molift - designed for life

At times with increasing demand for efficiency in the healthcare sector and a growing focus on patient needs, it is easy to neglect the needs of the carer. The need to give better help and caring avoiding heavy lifting situations, avoiding putting your own back at risk, without being completely exhausted after work and still have some surplus energy when the workday is over.

This is the insight that drives us at Molift. The company which was started 20 years ago by a man who through his own experience saw the need for lifting and moving patients in an effective and comfortable way – for both patient and carer.

We have since then delivered products that compliment the natural movement of the body, that is both simple and intuitive to use. Products designed to give the patient a better life in addition to providing the carer with real support and optimum work enjoyment.

This is the true meaning of; Molift – designed for life

molift

Molift Inc., 5008 W. Linebaugh Ave., Suite 60
Tampa, Florida 33624

Tel: (813) 969-2213, Toll free No.: 877-665-4387, www.molift.com